



**GENEVA SCIENTIFIC**  
LAB EQUIPMENT SOLUTIONS

**BEDDING DISPENSER - SBD**



Model	SBD
Exterior Dimensions (WxDxH)	39.37" x 29.53" x 76.77"
Noise Output Level (Measured as defined in ISO standard 3746)	Below 70 dBA (limit as per Decree 91/277)
Electrical	110V/60Hz 208V/60Hz 220V/60Hz

The Steelco semi automatic bedding dispenser SBD easy is a standalone unit dedicated to accurate dispensing of any kind of free flow bedding.

The unit requires no installation, it is plug-and-play with a total height limited to 76.77" (1950mm) to allow an easy access from a standard doorway or lift.

In order to protect the operator from allergene dust exposure, the unit integrates a high performance air suction system that surrounds the dispensing ports.

The bedding dispenser is characterized by a high hygienic level thanks to the full stainless steel AISI 304 construction. The surface of the control panel is made of one piece glass, a Steelco exclusive in the market, it does not have any recess where dust can collect and allows the easiest surface cleaning.

The bedding dispenser is controlled by an industrial grade PLC and the HMI is Soft Touch type with graphic color LCD display.

## BEDDING DISPENSER - SBD

### STANDARD FEATURES

#### **Bedding Loading/Refilling**

\*Bedding is loaded manually at ergonomic height (33.46" / 850mm) on to the loading hopper. The hinged supporting grid for the filling of the cages can be lifted up for bedding refilling and is maintained open by gas pistons.

\*The clean bedding falls into the hopper (37 gal) and is then conveyed via a screw conveyor in dispensing silo (55.5 gal) that is positioned in the upper part of the dispenser. The system ensures that a full bottom tank will transfer any 'free flowing bedding' to the top hopper in 1 minute approximately.

\*The bedding dispenser is suitable also for conveyed filling.

#### **Bedding Dispensing**

\*The dosing system uses capacitive sensors to detect the presence of the cages in the correct position to start dispensing and dispenses automatically without having to touch a trigger bar or even come in contact with the dispenser.

\*The standard configuration includes two delivery points that can operate either independently (when filled by a single mouth), or connected (when you fill the larger cages or you choose the continuous delivery from both mouths), this allows minimum waste from the upper reservoir.

\*The distribution of the bedding takes place by gravity, a timer sets the amount of bedding to be discharged.

\*A system of rotating blades in the dispenser grants the repeatability of the amount of bedding dispensed. That amount is not influenced by the volume of bedding present in the silo.

#### **Air Suction System**

\*In order to protect the operator from dust exposure, the dispenser is equipped with an efficient air suction system. Suction ports are concentrated around the dispensing mouths and work continuously to avoid allergens and dust shedding.

\*Easy access to the high level air filters for maintenance.

#### **System HMI**

\*Full glass control panel with soft touch button areas and graphic LCD display

#### **System Monitoring**

\*Hopper empty sensor

\*Over fill sensor - preventing over filling

#### **Microprocessor Control System**

\*Industrial grade PLC. The electronic programmable microprocessor allows the use of up to 40 dispensing programs.

\*Top tank level status

\*For each program, the operator can set the parameters of independent or connected mouth bedding dispensing and discharge time.

\*Automatic sleep mode: after 5 minutes of non-use the bedding dispenser will go to sleep turning off all running parts until reactivated.

#### **Flexibility**

\*The bedding dispensing unit SBD is on casters so it can be easily moved and placed wherever needed.

### SAFETY FEATURES

#### **Safety Shut-offs**

\*Automatic safety shut-off interlocks on the hinged supporting grid.

\*Automatic safety shut-off interlocks on fan filter.

### CONSTRUCTION

\*Entirely made of AISI 304 stainless steel

\*External panels made of AISI 304 stainless steel Scotch Brite finish Ra<40 µm and tempered HST glass

\*Insulation: high performance acoustic insulation guards

### OPTIONAL FEATURES

\***Seismic Tie Down:** Anchors bedding dispenser to floor

\***Anti-Clog Kit:** Anti-clog kit for compatibility with non free flow types of bedding

\***Multipoint Dispensing Configuration:** Bedding dispenser configured with more than 2 dispensing ports



**GENEVA SCIENTIFIC**  
LAB EQUIPMENT SOLUTIONS



Geneva Scientific LLC  
P.O. Box 408  
Fontana, WI 53125

1-877-436-3827  
Fax: 262-245-6678  
Sales@Geneva-Scientific.com

Specifications are subject to change without notice.  
rev. 10/15