



**GENEVA SCIENTIFIC**  
LAB EQUIPMENT SOLUTIONS

**DROSOPHILA CHAMBER - DR-41NL**



Model	DR-41NL
Exterior Dimensions (WxDxH)	41" x 33.6" x 77.2"
Temp. Range (w/ lights on)	2-44±0.5
Interior Space	37.2 cu.ft.
Total Shelving Floor Area	40.8 sq.ft.
Maximum Growing Height	8"
Light Intensity (6" from lamps)	No Light
Number of Tiers	6

**Applications**

- \*Designed especially for drosophila rearing
- \*Many other applications exist for this product

**Percival's IntellusUltra Controller**

\*Percival Scientific has built a reputation of providing flexible, customized options for research scientists around the world. We have taken the philosophy to the next level with our improved IntellusUltra Controller. Now choose from the levels of functionality that meet your needs.

**Lighting System**

- \*No lights are contained inside the chamber

**Cabinet Construction**

- \*26-gauge smooth galvanized white side walls and top reinforced with 24-gauge backer plates
- \*Interior floor constructed of 24-gauge #304-4 stainless steel
- \*24-gauge smooth white galvanized exterior
- \*NSF compliant seam design

**Cabinet Construction cont.**

- \*Overall wall thickness 2"
- \*Foamed-in-place non-CFC insulation
- \*One 1¼" diameter access port on right hand wall
- \*Chamber floor equipped with floor drain and hose assembly
- \*Contains casters assembly and adjustable leveling legs to compensate for floor unevenness in the lab

**Insulation**

- \*Woodless construction using CFC free insulation (overall wall thickness is 2", ample insulation for maintenance of stated temperature range)

**Doors**

- \*One door opening of 36.8" x 57.5" providing full access to chamber interior (magnetic gasket provides a tight seal to door frame)

**Interior Space**

- \*32.7 cu.ft. with work area of 40.8 sq.ft. provided on six tiers

**Shelving**

- \*Six tiers of white epoxy coated steel wire shelving (each shelf is 27"D x 36.3"W)

## DROSOPHILA CHAMBER - DR-41NL

### Shelving cont.

- \*Shelves are supported by shelf clips allowing 1/2" vertical adjustments
- \*Maximum clearance between shelves is 8" per tier with all six shelves installed (and the water pan humidifier as standard)

### Finish

- \*Interior and exterior painted with high reflective, environmentally friendly, high temperature baked white powder coating

### Refrigeration

- \*Self-contained air-cooled condensing unit with hot gas bypass system for continuous compressor operation, extended life and close temperature control (this continuous running condensing unit ensures precise temperature control by alternating cycling refrigerant and hot gas to coil; this also prolongs life of compressor, and eliminates risk of ice build up in coil)
- \*Top mounted refrigeration system allows easy access for maintenance (e.g. cleaning)
- \*As heat is rejected, it rises and is dissipated into the room without having any effect on inside temperature of cabinet
- \*Solenoid valves have extended stem for quiet and long life operation
- \*Ceiling mounted evaporator coils incorporate twin air circulation fans in aluminum housing (heat rejection to ambient with standard chamber: 3082 BTU/hr)
- \*Coil is phenolic coated to prolong the life of chamber from the acidic environment that can be created by drosophila

### Temperature Range

- \*2-44C (+/- 0.5C)

### Temperature Safety Limit Controls

- \*Experiment Protection: Adjustable high and low temperature controls, audible alarms, and visual indicators provided
- \*Controls shut down all power to chamber, activating alarms (when the temperature returns to the normal range the system will automatically reset)

### Humidity Control

- \*Additive humidity control of higher than ambient to 85% (+/-10%) lights on for set temperatures between 15° to 30°C
- \*Humidity control of higher than ambient to 90% (+/-10%) lights off for set temperatures between 15° to 30°C
- \*Humidity system requires a de-mineralized water supply line which terminates to a 1/2" MPT connector

### Electrical Requirements

- \*Consult Geneva Scientific for electrical requirements and amperage draw.

### Popular options

- \*Dehumidifier with Sensor
- \*IntellusUltra Connect
- \*Android-Based Touch Screen
- \*CO2 Enrichment Package
- \*Self-Contained Water-Cooled Condensing Unit
- \*Dry Alarm Contacts
- \*Convenience Receptacles



**GENEVA SCIENTIFIC**  
LAB EQUIPMENT SOLUTIONS



Geneva Scientific LLC  
P.O. Box 408  
Fontana, WI 53125

1-877-436-3827  
Fax: 262-607-6497  
Sales@Geneva-Scientific.com

Specifications are subject to change without notice.  
rev. 03/22