



# GENEVA SCIENTIFIC

LAB EQUIPMENT SOLUTIONS

## INCUBATOR - I-30L



Model	I-30L
Exterior Dimensions (WxDxH)	31" x 23.8" x 46.2"
Temp. Range (w/ lights on)	5-44±0.5
Interior Space	9.6 cu.ft.
Total Shelving Floor Area	9.1 sq.ft.
Maximum Growing Height	9.2"
Light Intensity (6" from lamps)	123 micromoles/m2/sec
Number of Tiers	2

### Applications

- \*This chamber is frequently used for Cyanobacteria, insects and simple/low cost seed germination
- \*Many other applications exist for this product

### Percival's IntellusUltra Controller

\*Percival Scientific has built a reputation of providing flexible, customized options for research scientists around the world. We have taken the philosophy to the next level with our improved IntellusUltra Controller. Now choose from the levels of functionality that meet your needs.

### Lighting System

- \*Linear LEDs horizontally mounted in pairs above each shelf
- \*Intensity programmable up to 125  $\mu\text{moles/m}^2/\text{s}$  of light irradiance measured @ 6" from LEDs
- \*Programming and control of the lighting is done via IntellusUltra real time controller
- \*Dimmable between 5-100% output

### Cabinet Construction

- \*Interior constructed of 22-gauge electro-zinc plated steel
- \*Exterior constructed of 18-gauge exterior electro-zinc plated steel
- \*Welded seams and joints on outer and inner shells
- \*Inner shell supported by non-compressing/non-thermal conducting material locking inner liner in place without a metal-to-metal bond to outer case
- \*Chamber is completely self-contained, suitable for stacking one above the other
- \*Overall wall thickness is 2"

### Insulation

- \*Woodless construction using CFC free insulation (overall wall thickness is 2", ample insulation for maintenance of stated temperature range)

### Door

- \*One door opening of 26.8" x 29.6" providing full access to chamber interior (magnetic gasket provides a tight seal to door frame)

### Interior Space

- \*9.6 ft<sup>3</sup> (0.3 m<sup>3</sup>) with a work area 9.1 ft<sup>2</sup> (0.8 m<sup>2</sup>) provided on two tiers

## INCUBATOR - I-30L

### Shelving

\*Two tiers of white epoxy coated steel wire shelving (each shelf is 16.5"D x 26.5"W)

\*Shelves are supported by shelf clips allowing ½" vertical adjustments

\*Maximum clearance between shelves is 9.2" per tier with both shelves installed

### Finish

\*Interior and exterior painted with high reflective, environmentally friendly, high temperature baked white powder coating

### Refrigeration

\*Self-contained air-cooled condensing unit with hot gas bypass system for continuous compressor operation, extended life and close temperature control (this continuous running condensing unit ensures precise temperature control by alternating cycling refrigerant and hot gas to coil; this also prolongs life of compressor, and eliminates risk of ice build up in coil)

\*Solenoid valves have extended stem for quiet and long life operation

\*Rear chamber wall mounted evaporator coil incorporates an air circulation fan (heat rejection to ambient with standard chamber: 1442 BTU/hr)

### Temperature Range

\* 2°-44°C (±0.5°C) lights off and 5°-44°C (±0.5°C) lights on

### Temperature Safety Limit Controls

\*Experiment Protection: Adjustable high and low temperature controls, audible alarms, and visual indicators provided

\*Controls shut down all power to chamber, activating alarms (when the temperature returns to the normal range the system will automatically reset)

### Electrical Requirements

\*Consult Geneva Scientific for electrical requirements and amperage draw.

### Options

\*Phenolic Coated Coils

\*Access Port with cover

\*Additive Humidity Control with Sensor

\*Dehumidifier with Sensor

\*IntellusUltra Connect

\*Android-Based Touch Screen

\*CO2 Enrichment Package

\*Self-Contained Water-Cooled Condensing Unit

\*Dry Alarm Contacts

\*Open loop dimmable lighting per tier

\*Convenience receptacles

\*Extended temperature ranges



**GENEVA SCIENTIFIC**  
LAB EQUIPMENT SOLUTIONS

